



Made of 18/8 Stainless Steel, Solid, Set of 10



All The Above Impression Trays are Available in These 5 Sizes




UPPER		U1	U2	U3	U4	U5
	A	62	60	53	54	52
B	83	75	73	69	68	


All Dimension are in mm

LOWER		L1	L2	L3	L4	L5
	A	63	61	57	60	53
B	82	81	78	74	71	


All Dimension are in mm




Perforated, Set of 10



Made of 18/8 Stainless Steel, Solid Full Denture Set of 8




All The Above Impression Trays are Available in These 4 Sizes




UPPER		XL	L	M	S
	A	70	67	64	60
B	79	76	70	73	

All Dimension are in mm




LOWER		XL	L	M	S
	A	64	64	60	57
B	86	70	73	70	

All Dimension are in mm




Perforated Full Denture Set of 8

Made of 18/8 Stainless Steel, Solid Set of 6



All The Above Impression Trays are Available in These 4 Sizes





UPPER		U1	U2	U3
	A	67	64	57
B	81	71	70	

All Dimension are in mm

LOWER		L1	L2	L3
	A	60	57	54
B	81	79	76	

All Dimension are in mm




Perforated Set of 6

Made of 18/8 Stainless Steel, English Pattern Set of 8



All The Above Impression Trays are Available in These 4 Sizes




<b>UPPER</b>	Size	A	B	C
	0	67	50	12
	1	72	51	13
	2	78	55	15
	3	80	58	15


All Dimension are in mm

<b>LOWER</b>	Size	A	B	C
	0	70	50	14
	1	77	51	15
	2	79	55	15
	3	80	60	17


All Dimension are in mm



Made of 18/8 Stainless Steel, Solid Full Denture Set of 8



All The Above Impression Trays are Available in These 4 Sizes





UPPER	XL	L	M	S
A	70	67	64	60
B	79	76	70	73

All Dimension are in mm


LOWER	XL	L	M	S
A	64	64	60	57
B	86	70	73	70

All Dimension are in mm



Made of 18/8 Stainless Steel, Perforated Full Denture Set of 8

Made of 18/8 Stainless Steel, Solid Full Denture Set of 8




Upper Sizes: 0,1,2,3

Depth of Tray: 15,9 mm, Height of Plate: 12,5 mm


Lower Sizes: 0,1,2,3

Depth of Tray: 12,9 mm



Made of 18/8 Stainless Steel, Perforated Full Denture Set of 8

Impression Trays For Crown & Bridge Work



**PERFORATED**

**SOLID**



2931-A  
For upper right jaw  
For lower left jaw

2931-B  
For upper left jaw  
For lower right jaw



2931-C  
For upper right jaw  
For lower left jaw

2931-D  
For upper left jaw  
For lower right jaw

**Available in 4 sizes**

- 1.....44 mm
- 2.....55 mm
- 3.....60 mm
- 4.....65 mm

**Available in 4 sizes**

- 1.....44 mm
- 2.....55 mm
- 3.....60 mm
- 4.....65 mm