

Plaster Instruments

Plaster Knives



with ebony handle

07-9001
HOPKINS
200 mm



07-9002
ESMARCH
180 mm

Plaster Saw



07-9003
ENGEL
160 mm

Bandage and First-Aid scissors



07-9004 - 07-9012 145mm

07-9004 black	07-9008 yellow
07-9005 red	07-9009 white
07-9006 blue	07-9010 orange
07-9007 light green	07-9011 violet
	07-9012 bright pink



07-9013 - 07-9021 180mm

07-9013 black	07-9017 yellow
07-9014 red	07-9018 white
07-9015 blue	07-9019 orange
07-9016 light green	07-9020 violet
	07-9021 bright pink

Plaster Scissors


For small plaster dressings



Plaster Shears

Special features:

The construction of the open jaws avoids tiresome clogging by plaster remnants and allows furthermore an easy cleaning. The shears do not glide off the plaster owing to the deep teething on the jaws.
The knives with hollow grinding have a better cutting ability and therefore less force is required.
The knives are guided by ball bearings.



07-9029
STILLE
260 mm



07-9030
STILLE
360 mm

STILLE Plaster Shears



Blade with reinforced back



1/3

07-9032
STILLE
230 mm

07-9033
Spare blade



1/3

07-9034
STILLE
370 mm

07-9035
Spare blade



1/3

07-9036
STILLE
240 mm

07-9037
Spare blade

Plaster Spreader, Plaster Breakers



1/3

07-9038
HENNIG
270 mm




1/3

07-9039
MOD.USA
230 mm



1/3

07-9040
290 mm



with removable jaws



1/3

07-9041 - 07-9042
WOLFF
07-9041 180 mm
07-9042 240 mm
plaster breakers