

# Vacuum Currettes for Aspiration Lipectomy, Suction Instruments



**07-9043 - 07-9048**  
BIEMER  
07-9043 ∅ 4mm  
07-9044 ∅ 6mm  
07-9045 ∅ 8mm  
07-9046 ∅ 10mm  
07-9047 ∅ 12mm  
07-9048 ∅ 14mm  
Currettes for aspiration lipectomy, graduated, slightly curved

**07-9049 - 07-9050**  
BIEMER  
07-9049 ∅ 4 mm  
07-9050 ∅ 6 mm  
Currettes for aspiration lipectomy, with three apertures, slightly curved malleable

**07-9051 - 07-9053**  
BIEMER  
07-9051 ∅ 8 mm  
07-9052 ∅ 10 mm  
07-9053 ∅ 12 mm  
Currettes for aspiration lipectomy, with three apertures, slightly curved malleable


**07-9054**  
For aspirating fluid from the excretory ducts of large glands and for keeping the operating area clear  
170 mm

**07-9055**  
Trocar  
For puncturing  
150 mm

**07-9056**  
For clearing large quantities of fluid  
190 mm


# Suction Instruments

**07-9057** 150mm  
Suction handle, with vacuum regulator and finger cut-off, brass, mirror-finish chromium-plated




olive  
Ø 9 mm


**07-9058** 170mm  
Suction handle, with hand regulation valve and quick-shutter, brass, mirror-finish chromium-plated



★ Both suction handles have bayonet fitting olive and can be used for the suction cannulae. Ø 9 mm




**07-9059**  
Trocar  
For puncturing  
170 mm



**07-9060**  
Hollow needle,  
sharp  
150 mm



**07-9061**  
double curved  
For general  
surgery  
230 mm



**07-9062**  
For tonsillectomy,  
oral and laryngeal  
operations  
190 mm




Ø 11 mm  
Ø 7,5 mm

**07-9063**  
For tonsillectomy,  
oral and laryngeal  
operations  
220 mm

★ With bayonet fitting to be used with suction handles (07-9057 - 07-9058)

# Suction Instruments



**07-9064**  
For tonsillectomy,  
oral-and laryngeal  
operations  
180 mm

**07-9065**  
Suction cannula  
for tonsillectomy,  
oral-and laryngeal  
operations  
220 mm

**07-9066**  
Suction cannula  
for tonsillectomy,  
oral-and laryngeal  
operations  
180 mm

**07-9067**  
S-shaped blunt cannula  
for aspirating fluid from  
the excretory ducts of  
large glands and for  
keeping the operating  
area clear  
170 mm


**07-9068**  
POOL  
Suction cannula,  
with outer tube  
to be screwed of  
225 mm

**07-9069**  
POOL  
Suction cannula,  
with outer tube  
to be screwed off  
225 mm

# Suction Instruments




# Suction Instruments



07-9076 Ø 1,5 mm  
07-9077 Ø 2,0 mm  
07-9078 Ø 2,5 mm  
07-9079 Ø 3,0 mm

**07-9076 - 07-9079**  
**PLESTER**  
with finger cut-off and  
stylet, thumb piece  
parallel to curve, with  
LUER hub  
195 mm




07-9080 Ø 2,0 mm  
07-9081 Ø 2,5 mm  
07-9082 Ø 3,0 mm  
07-9083 Ø 4,0 mm  
07-9084 Ø 5,0 mm


**07-9080 - 07-9084**  
**FERGUSSON**  
with finger cut-off and  
stylet, thumb piece  
90° to curve, with  
LUER hub  
180 mm



**07-9085** Ø 4mm  
**FERGUSSON**  
insulated, with finger cut-off  
and stylet, thumb piece 90°  
to curve, with LUER hub  
180 mm




07-9086 Ø 3mm  
07-9087 Ø 5mm  
**FRAZIER**  
with finger  
cut-off and stylet, thumb  
piece 90° to curve,  
with LUER hub  
180 mm



**07-9088**  
5,0 mm diam., Tip  
6,5 mm wide, Nasal  
suction tube, with  
LUER hub  
180 mm




07-9089 4,5 mm  
07-9090 5,5 mm  
Nasal suction  
tubes  
175 mm



**07-9091**  
**GORNEY**  
Suction dissector  
180 mm



07-9092 4 mm wide  
07-9093 5 mm wide  
Suction raspatories,  
sharp, bayonet-shaped,  
with stylet  
200 mm



**07-9094**  
**LOPEZ-REINKE**  
Tonsil dissector with  
suction device  
235 mm

# Suction Instruments

**07-9095**  
WULLSTEIN  
Suction tube, conical, silver, olive  
with LUER hub  
160 mm

**07-9096** 1,5 mm  
**07-9097** 2,0 mm  
WULLSTEIN  
Suction tube, cylindrical, silver, with  
stylet olive with LUER hub  
160 mm

olive Ø8 mm

olive Ø8 mm

Usable length

YASARGIL usable length diam.		YASARGIL usable length diam.	
<b>07-9098</b>	130 mm 1,5 mm	<b>07-9106</b>	130 mm 2,5 mm
<b>07-9099</b>	150 mm 1,5 mm	<b>07-9107</b>	150 mm 2,5 mm
<b>07-9100</b>	180 mm 1,5 mm	<b>07-9108</b>	180 mm 2,5 mm
<b>07-9101</b>	220 mm 1,5 mm	<b>07-9109</b>	220 mm 2,5 mm
<b>07-9102</b>	130 mm 2,0 mm	<b>07-9110</b>	130 mm 3,5 mm
<b>07-9103</b>	150 mm 2,0 mm	<b>07-9111</b>	150 mm 3,5 mm
<b>07-9104</b>	180 mm 2,0 mm	<b>07-9112</b>	180 mm 3,5 mm
<b>07-9105</b>	220 mm 2,0 mm	<b>07-9113</b>	220 mm 3,5 mm
		<b>07-9114</b>	180 mm 4,5 mm

olive  
Ø6,5-9 mm  
with LUER hub

**07-9115 - 07-9117**  
LANDOLT  
150 MM  
diam Ø 3 mm

**07-9115** 8 mm / curved  
**07-9116** 11mm / curved  
**07-9117** 15mm / curved

**07-9118** 500 mm  
**07-9119** 350 mm  
with rubber tip,  
with LUER hub

olive Ø9 mm

**07-9120** 2,0 mm diam  
**07-9121** 2,5 mm diam  
**07-9122** 3,0 mm diam  
with finger cut-off,  
thumb piece 90° to  
curve, with LUER hub  
250 mm

olive Ø6-9 mm

# Combined Suction and Irrigation Cannulae, Coagulation Suction Cannulae

